Vitamin E

Vitamin E is a fat soluble vitamin of which α tocopherol has the highest biological activity. It is found in vegetables, seed oils and nuts. It is absorbed in the small intestine with about half passing to the liver and half passing into the lymph. Within the blood it is transported bound to LDL and HDL. Vitamin E accumulates in the liver, muscle and adipose tissue.

Vitamin E is an important chain-breaking antioxidant in cell membranes. It also acts to modulate signal transduction pathways.

Deficiency

Deficiency is rare but may occur in

- Premature and newborn babies
- Conditions associated with fat malabsorption
- Patients on parenteral nutrition
- Abetalipoproteinaemia
- Homozygous hypobetalipoproteinaemia
- Familial isolated Vitamin E deficiency

This may result in *peripheral neuropathy* and *spinocerebellar syndrome*. Note that newborn babies have low serum Vitamin E concentration compared with adults, but may not have true Vitamin E deficiency ie with symptoms.

Toxicity

Relative to other fat-soluble vitamins, Vitamin E is safe. Few side effects have been reported even on very high doses.

Reference ranges

0 - <1yr: 5-50 μmol/L Source: CALIPER Clin Biochem (2014) 47(9):812-815

1-6 yrs:7-21 μmol/L7-12 yrs:10-21 μmol/L13-18 yrs:13-24 μmol/LSource: Clin Chem 1988 34(8) 1625-1628

Adult: $11.6 - 35.5 \mu mol/L$ Source: Local laboratory data performed as part of FIMLS dissertation 1989

Vitamin E is bound to lipoproteins and the results of Vitamin E analysis are reported as a molar ratio with serum lipids if elevated (and if enough sample supplied)

Ratio = Vitamin E (µmol/L) / (cholesterol + triglyceride) (mmol/L) The ratio is a good indicator of tissue stores

Paediatric Vitamin E/Lipid ratioNo molar ratio available for < 1 year</td>1-6 yrs3-57-12 yrs2-513-19 yrs2-4Source: Clin Chem 1988 34(8) 1625-1628

Adult Vitamin E/Lipid ratio 3.9-5.9 Source: Local Study data

Specimen type

Serum (also lithium heparin or EDTA plasma) Protect from light. Minimum volume 50 μ L – done on dilution (ideally > 100 μ L) NOT haemolysed (haemolysed samples can falsely lower results)

Storage Separate serum or plasma asap and freeze

Transport

First class post, ambient temperature. Protect from light

Address for specimens

Department of Clinical Biochemistry Rotherham Hospital Moorgate Road Rotherham, S60 2UD

Cost

Contact - annettedavis-green@nhs.net

Turnaround Weekly HPLC assay

Accreditation

Accredited to UKAS ISO15189

External QA

UKNEQAS

Contact persons

Consultant Clinical Scientist rgh-tr.biochemistry@nhs.net Tel 01709 820000 (Hospital) Tel 01709 424051 (Secretary)

(September 2021 V5)

